




Authorization chain – every gate must pass, each independently revocable

- 1 · Customer consent
- 2 · Product license
- 3 · Engineer MFA
- 4 · RBAC role
- 5 · Class-II review
- 6 · PHO review

Deny-by-default

The host carries 0 standing access policies
No session ⇒ no policy ⇒ no certificate ⇒ no access. Close ⇒ revoked ≤ 1 heartbeat.
The console is the gate.

Why "Authenticated using none" is not a backdoor

That `none` is only the client→WARP-proxy leg. Authentication to the host is a short-lived certificate signed by the Cloudflare SSH CA, validated by the host's TrustedUserCAKeys – no standing secret on the host, and the cert exists only under an authorized session.

Provisioning

The session is recorded and tied to a service ticket.

- no standing access
- session-bound
- after close, revoked